

101896GR - 8.1

QUANSYS
 BIOSCIENCES
CERTIFICATE OF ANALYSIS


Quansys Biosciences
 365 N 600 W
 Logan, UT 84321
 (435) 752-0531

Q-Plex™ Human Pneumococcal IgG (9-Plex)

REF	Part Number: 482149HU	Software Product Code: HPNM231211-9.3 
LOT	Lot Number: HPNK231211 21925 001	
	Expiration Date: 2026-JAN-07	

Reconstitute Calibrator using 500 µL of Sample Diluent

Component	Part Number	Lot Number	Recon. Volume	Exp. Date
Plate	482152HU	HPNP250128	N/A	2026-JAN-28
Sample Diluent	481157HU	HPNS250107	N/A	2026-JAN-07
Calibrator	481155HU	HPNM231211	500 µL	2026-MAR-07
Detection	482156HU	HPND250121	N/A	2026-JAN-21
Substrate A	101159GR	QNSA250113	N/A	2026-DEC-09
Substrate B+	101120GR	QNSB+250113	N/A	2026-DEC-09
Wash Buffer	101158GR	QNSW250106	N/A	2027-JAN-06

**Q-Plex™** ARRAYS

Spot #	Analyte	Cal Conc	Units	LLD	LLOQ	ULOQ
1	CWPS	N/A	N/A	N/A	N/A	N/A
2	2	2.25	µg/mL	0.0011	0.0022	0.56
3	10A	1.19	µg/mL	0.00058	0.0012	0.074
4	11A	0.47	µg/mL	0.00023	0.00046	0.12
5	15B	1.55	µg/mL	0.00076	0.0015	0.39
6	17F	0.73	µg/mL	0.0014	0.0028	0.18
7	19A	1.28	µg/mL	0.0025	0.0050	0.32
8	20	0.92	µg/mL	0.00045	0.00090	0.23
9	22F	0.87	µg/mL	0.0017	0.0034	0.22
10	33F	0.94	µg/mL	0.00046	0.00092	0.059
11	Control Spot	N/A	N/A	N/A	N/A	N/A
12	REF	N/A	N/A	N/A	N/A	N/A

QC Release:

Date:

01 MAR 25

QA Release:

for Research Use Only

Date:

01 MAR 25