

101896GR - 8.1

QUANSYS BIOSCIENCES
CERTIFICATE OF ANALYSIS


Quansys Biosciences
 365 N 600 W
 Logan, UT 84321
 (435) 752-0531

Q-Plex™ Human T Cell Cytokine HS (9-Plex)

REF	Part Number: 111933HU	Software Product Code: HCHM200803-TC9 
LOT	Lot Number: HTCK200803 22047 001	
	Expiration Date: 2026-MAR-19	

Reconstitute Calibrator using 750 µL of Sample Diluent

Component	Part Number	Lot Number	Recon. Volume	Exp. Date
Plate	111952HU	HTCP240619	N/A	2026-JUN-19
Sample Diluent	101104HU	HQ6S240805	N/A	2026-AUG-05
Calibrator	112455HU	HCHM200803	750 µL	2027-MAR-18
Detection	112456HU 01	HCHDQ6250214	6 mL	2027-FEB-14
Streptavidin HRP	101173GR	QNSN241105	N/A	2026-MAY-05
Substrate A	101159GR	QNSA240906	N/A	2026-MAR-19
Substrate B+	101120GR	QNSB+240906	N/A	2026-MAR-19
Wash Buffer	101158GR	QNSW250106	N/A	2027-JAN-06

**Q-Plex™ ARRAYS**

Spot #	Analyte	Cal Conc	Units	LLD	LLOQ	ULOQ
1	IFN γ	93	pg/mL	0.064	0.13	93
2	IL-2	47	pg/mL	0.032	0.064	47
3	IL-4	18	pg/mL	0.012	0.025	18
4	IL-5	260	pg/mL	0.18	0.36	260
5	IL-10	2320	pg/mL	1.59	3.18	2320
6	IL-13	290	pg/mL	0.20	0.40	290
7	IL-17	4500	pg/mL	3.09	6.17	4500
8	IL-21	260	pg/mL	0.18	0.36	260
9	TNF α	640	pg/mL	0.44	0.88	640
10	REF	N/A	N/A	N/A	N/A	N/A

QC Release:

R. Wa
m. Ay

Date:

21MAR25

QA Release:

Date:

21MAR25

for Research Use Only