

101896GR - 8.1

QUANSYS
 BIOSCIENCES
CERTIFICATE OF ANALYSIS


Quansys Biosciences
 365 N 600 W
 Logan, UT 84321
 (435) 752-0531

Q-Plex™ Porcine Cytokine HS (4-Plex)

REF	Part Number: 119149PC	Software Product Code: PCSM170306CY 
LOT	Lot Number: PCSK170306 22004 001	
	Expiration Date: 2026-MAR-19	

Reconstitute Calibrator using 500 µL of Sample Diluent


Component	Part Number	Lot Number	Recon. Volume	Exp. Date
Plate	119152PC	PCSP250128	N/A	2027-JAN-28
Sample Diluent	119157PC	PCSS250103	N/A	2027-JAN-03
Calibrator	119155PC	PCSM170306	500 µL	2027-MAR-10
Detection	119156PC	PCSD250107	N/A	2027-JAN-07
Streptavidin HRP	101173GR	QNSN241105	N/A	2026-MAY-05
Substrate A	101159GR	QNSA241108	N/A	2026-MAR-19
Substrate B+	101120GR	QNSB+241108	N/A	2026-MAR-19
Wash Buffer	101158GR	QNSW250106	N/A	2027-JAN-06



Spot #	Analyte	Cal Conc	Units	LLD	LLOQ	ULOQ
1	IL-1β	1600	pg/mL	N/A	N/A	N/A
2	IL-6	2000	pg/mL	N/A	N/A	N/A
3	IL-8	2300	pg/mL	N/A	N/A	N/A
4	TNFα	1350	pg/mL	N/A	N/A	N/A
5	REF	N/A	N/A	N/A	N/A	N/A

Q-Plex™ ARRAYS


QC Release:



Date:

12MAR25

QA Release:



Date:

12MAR25

for Research Use Only