

101896GR - 8.1

QUANSYS BIOSCIENCES
CERTIFICATE OF ANALYSIS


Quansys Biosciences
 365 N 600 W
 Logan, UT 84321
 (435) 752 - 0531

Q-Plex™ NHP Custom 6-Plex

REF	Part Number: 107749GR-06	Software Product Code: NHPM220606-SB6 
LOT	Lot Number: NCUK220606 22116 001	
	Expiration Date: 2026-MAR-19	

Reconstitute Calibrator 1 using 500 µL of Complete Sample Diluent
Reconstitute Calibrator 2 using 350 µL of Calibrator 1

Component	Part Number	Lot Number	Recon. Volume	Exp. Date
Plate	101852GR	NCUP250328	N/A	2026-SEP-28
Sample Diluent	101104HU	HQ6S240805	N/A	2026-AUG-05
Competitor	101890GR	NCUC250327	10 mL	2026-SEP-30
Calibrator 1	109255MB	NHPM220606	500 µL	2026-SEP-30
Calibrator 2	101855GR	NCUM250326	350 µL	2026-SEP-30
Detection	101856GR	NCUD250325	6 mL	2026-SEP-25
Streptavidin HRP	101173GR	QNSN241105	N/A	2026-MAY-05
Substrate A	101159GR	QNSA241108	N/A	2026-MAR-19
Substrate B+	101120GR	QNSB+241108	N/A	2026-MAR-19
Wash Buffer	101158GR	QNSW250106	N/A	2027-JAN-06



Spot #	Analyte	Cal Conc	Units	LLD	LLOQ	ULOQ
1	*CRP	3.90	mg/L	0.0080	0.016	3.90
2	IFN γ	820	pg/mL	0.56	1.12	820
3	IL-1 β	3660	pg/mL	2.51	5.02	3660
4	IL-6	910	pg/mL	0.62	1.25	910
5	IL-10	2460	pg/mL	1.69	3.37	2460
6	TNF α	190	pg/mL	0.13	0.26	190
7	REF	N/A	N/A	N/A	N/A	N/A

Q-Plex™ ARRAYS

*Competitive Assay

QC Release:

Date: 31 MAR 25

QA Release:

Date: 31 MAR 25

for Research Use Only