

101896GR - 8.1

QUANSYS BIOSCIENCES
CERTIFICATE OF ANALYSIS


Quansys Biosciences
365 N 600 W
Logan, UT 84321
(435) 752 - 0531

Q-Plex™ Human Custom 2-Plex

REF	Part Number: 107749GR-02	Software Product Code: HCUM200203-2P 
LOT	Lot Number: HCUK180417 22363 001	
	Expiration Date: 2026-JUL-24	

Reconstitute 2 vials of Calibrator 1 using 200 µL of Sample Diluent**Reconstitute Calibrator 2 using 200 µL of combined Calibrator 1 vials**


Component	Part Number	Lot Number	Recon. Volume	Exp. Date
Plate	101852GR	HCUP250520	N/A	2027-MAY-20
Sample Diluent	101104HU	HQ6S240805	N/A	2026-AUG-05
Assay Diluent	462183HU	HADV240805	N/A	2026-AUG-05
Competitor 1	462190HU	HCUC211007	10 mL	2027-MAY-29
Competitor 2	101890GR	HCUC180417	10 mL	2027-MAY-29
Calibrator 1	462155HU	HCUM200203	200 µL	2027-MAY-29
Calibrator 2	101855GR	HCUM191007	200 µL	2027-MAY-29
Streptavidin HRP	101173GR	QNSN250124	N/A	2026-JUL-24
Substrate A	101159GR	QNSA250113	N/A	2026-DEC-09
Substrate B+	101120GR	QNSB+250113	N/A	2026-DEC-09
Wash Buffer	101158GR	QNSW250106	N/A	2027-JAN-06



Spot #	Analyte	Cal Conc	Units	LLD	LLOQ	ULOQ
1	Adiponectin	25000	ng/mL	17.15	34.29	25000
2	VCAM-1	1213700	pg/mL	2497.33	4994.65	1213700
3	REF	N/A	N/A	N/A	N/A	N/A

Q-Plex™ ARRAYS

QC Release:



Date:

29 MAY 25

QA Release:



Date:

29 MAY 25

for Research Use Only