

101896GR - 8.1

QUANSYS
BIOSCIENCES
CERTIFICATE OF ANALYSIS


Quansys Biosciences
365 N 600 W
Logan, UT 84321
(435) 752-0531

Q-Plex™ Human Custom 5-Plex

REF	Part Number: 107749GR-05	Software Product Code:
LOT	Lot Number: RESK250528 22419 001	HRBM250528 
	Expiration Date: 2026-MAY-21	

Reconstitute Calibrator using 300 µL of Sample Diluent

Component	Part Number	Lot Number	Recon. Volume	Exp. Date
Plate	101152GR	RESP250521	N/A	2026-MAY-21
Sample Diluent	101104HU	HQ6S240805	N/A	2026-AUG-05
Calibrator	101855GR	HRBM250528	300 µL	2026-JUN-09
Control 1	101131GR	HRBQ250528-1	500 µL	2026-JUN-09
Control 2	101132GR	HRBQ250528-2	500 µL	2026-JUN-09
Detection	101856GR	HRBD250604	N/A	2026-JUN-04
Streptavidin HRP	101173GR	QNSN250124	N/A	2026-JUL-24
Substrate A	101159GR	QNSA250227	N/A	2026-DEC-09
Substrate B+	101120GR	QNSB+250227	N/A	2026-DEC-09
Wash Buffer	101158GR	QNSW250113	N/A	2027-JAN-13



Spot #	Analyte	Cal Conc	Units
1	PKM2	357	Units/mL
2	TIMP-1	253	ng/mL
3	IGFBP2	842	ng/mL
4	DKK3	32	ng/mL
5	BDNF	29	ng/mL
6	REF	N/A	N/A

Q-Plex™ ARRAYS

QC Release:

Date: 10 JUN 25

QA Release:

Date: 10 JUN 25*for Research Use Only*