

101896GR - 8.1

QUANSYS
 BIOSCIENCES
CERTIFICATE OF ANALYSIS


Quansys Biosciences
 365 N 600 W
 Logan, UT 84321
 (435) 752 - 0531

Q-Plex™ Human Custom 3-Plex

REF	Part Number: 107749GR-03	Software Product Code: HA1M220719-3P 
LOT	Lot Number: HCUK220719 22529 001	
	Expiration Date: 2026-MAR-19	

Reconstitute Calibrator using 650 µL of Sample Diluent


Component	Part Number	Lot Number	Recon. Volume	Exp. Date
Plate	101852GR	HCUP250721	N/A	2027-JUL-21
Sample Diluent	101104HU	HQ6S240805	N/A	2026-AUG-05
Calibrator	102655HU	HA1M220719	650 µL	2027-JUL-25
Detection	119056HU 01	HCSDQ6250516	N/A	2027-MAY-16
Streptavidin HRP	101173GR	QNSN241105	N/A	2026-MAY-05
Substrate A	101159GR	QNSA240906	N/A	2026-MAR-19
Substrate B+	101120GR	QNSB+240906	N/A	2026-MAR-19
Wash Buffer	101158GR	QNSW241113	N/A	2026-NOV-13



Spot #	Analyte	Cal Conc	Units	LLD	LLOQ	ULOQ
1	IL-1β	11000	pg/mL	7.54	15.09	11000
2	IL-10	7870	pg/mL	5.40	10.80	7870
3	TNFα	5490	pg/mL	3.77	7.53	5490
4	REF	N/A	N/A	N/A	N/A	N/A

Q-Plex™ ARRAYS

QC Release:



Date:

25JUL25

QA Release:



Date:

25JUL25

for Research Use Only