

101896GR - 8.1

QUANSYS
 BIOSCIENCES
CERTIFICATE OF ANALYSIS


Quansys Biosciences
 365 N 600 W
 Logan, UT 84321
 (435) 752-0531

Q-Plex™ Human Angiogenesis (9-Plex)

REF	Part Number: 150233HU	Software Product Code: HA1M220719-AG 
LOT	Lot Number: HAGK220719 22724 001	
	Expiration Date: 2026-JAN-05	

Reconstitute Calibrator using 650 µL of Sample Diluent

Component	Part Number	Lot Number	Recon. Volume	Exp. Date
Plate	150252HU	HAGP250818	N/A	2027-AUG-18
Sample Diluent	101104HU	HQ6S240805	N/A	2026-AUG-05
Calibrator	102655HU	HA1M220719	650 µL	2027-SEP-17
Detection	150256HU	HAGDQ6240105	N/A	2026-JAN-05
Streptavidin HRP	101173GR	QNSN241105	N/A	2026-MAY-05
Substrate A	101159GR	QNSA240531	N/A	2026-MAR-19
Substrate B+	101120GR	QNSB+240531	N/A	2026-MAR-19
Wash Buffer	101158GR	QNSW241030	N/A	2026-OCT-30

**Q-Plex™** ARRAYS

Spot #	Analyte	Cal Conc	Units	LLD	LLOQ	ULOQ
1	ANG-2	9070	pg/mL	6.22	12.44	9070
2	FGF	4530	pg/mL	3.11	6.21	4530
3	HGF	13300	pg/mL	9.12	18.24	13300
4	IL-8	1880	pg/mL	1.29	2.58	1880
5	PDGF	3910	pg/mL	2.68	5.36	3910
6	TIMP-1	12300	pg/mL	25.31	50.62	12300
7	TIMP-2	18100	pg/mL	12.41	24.83	18100
8	TNFα	5490	pg/mL	3.77	7.53	5490
9	VEGF	1770	pg/mL	1.21	2.43	1770
10	REF	N/A	N/A	N/A	N/A	N/A

QC Release:



Date:

25SEP25

QA Release:



Date:

25SEP25*for Research Use Only*