

101896GR - 8.1

**QUANSYS**   
BIOSCIENCES  
**CERTIFICATE OF ANALYSIS**


Quansys Biosciences  
365 N 600 W  
Logan, UT 84321  
(435) 752-0531

**Q-Plex™ Human Custom Singleplex**

|                                                                                   |                                  |                                                                                                                                        |
|-----------------------------------------------------------------------------------|----------------------------------|----------------------------------------------------------------------------------------------------------------------------------------|
| <b>REF</b>                                                                        | Part Number: 107749GR-01         | Software Product Code:<br><b>FCHM190715-FSH</b><br> |
| <b>LOT</b>                                                                        | Lot Number: HCUK190715 22697 001 |                                                                                                                                        |
|  | Expiration Date: 2026-AUG-21     |                                                                                                                                        |


**Reconstitute Calibrator using 550 µL of Sample Diluent**

| Component        | Part Number | Lot Number  | Recon. Volume | Exp. Date   |
|------------------|-------------|-------------|---------------|-------------|
| Plate            | 101852GR    | HCUP250821  | N/A           | 2026-AUG-21 |
| Sample Diluent   | 101122HU    | HQ16250825  | N/A           | 2026-AUG-25 |
| Calibrator       | 332155HU    | FCHM190715  | 550 µL        | 2027-APR-17 |
| Detection        | 101856GR    | HCUD250828  | N/A           | 2026-AUG-28 |
| Streptavidin HRP | 101173GR    | QNSN250610  | N/A           | 2026-DEC-10 |
| Substrate A      | 101159GR    | QNSA250508  | N/A           | 2026-DEC-09 |
| Substrate B+     | 101120GR    | QNSB+250508 | N/A           | 2026-DEC-09 |
| Wash Buffer      | 101158GR    | QNSW250113  | N/A           | 2027-JAN-13 |



| Spot # | Analyte | Cal Conc | Units  | LLD  | LLOQ | ULOQ |
|--------|---------|----------|--------|------|------|------|
| 1      | FSH     | 193      | mIU/mL | 0.40 | 0.79 | 193  |
| 2      | REF     | N/A      | N/A    | N/A  | N/A  | N/A  |

**Q-Plex™** ARRAYS

QC Release: Date: 20 OCT 25QA Release: Date: 20 OCT 25**for Research Use Only**