

101896GR - 8.1

QUANSYS
BIOSCIENCES
CERTIFICATE OF ANALYSIS


Quansys Biosciences
365 N 600 W
Logan, UT 84321
(435) 752-0531

Q-Plex™ Human Custom 2-Plex

REF	Part Number: 107749GR-02	Software Product Code: HA1M220719-2P 
LOT	Lot Number: HCUK220719 22805 001	
	Expiration Date: 2026-MAR-19	

Reconstitute Calibrator 1 using 650 µL of Sample Diluent
Reconstitute Calibrator 2 using 400 µL of Calibrator 1

Component	Part Number	Lot Number	Recon. Volume	Exp. Date
Plate	101852GR	HCUP251002	N/A	2027-OCT-02
Sample Diluent	101104HU	HQ6S240805	N/A	2026-AUG-05
Calibrator 1	102655HU	HA1M220719	650 µL	2027-OCT-13
Calibrator 2	102855HU	HA2M200303	400 µL	2027-OCT-13
Detection	101856GR	HCUD251007	N/A	2027-OCT-07
Streptavidin HRP	101173GR	QNSN241105	N/A	2026-MAY-05
Substrate A	101159GR	QNSA240531	N/A	2026-MAR-19
Substrate B+	101120GR	QNSB+240531	N/A	2026-MAR-19
Wash Buffer	101158GR	QNSW250113	N/A	2027-JAN-13



Spot #	Analyte	Cal Conc	Units	LLD	LLOQ	ULOQ
1	MIG	8860	pg/mL	6.08	12.15	8860
2	IP-10	3240	pg/mL	2.22	4.44	1080
3	REF	N/A	N/A	N/A	N/A	N/A

Q-Plex™ ARRAYS

QC Release:



Date: 13 OCT 25

QA Release:



Date: 13 OCT 25

for Research Use Only