

101896GR - 8.1

QUANSYS
BIOSCIENCES
CERTIFICATE OF ANALYSIS


Quansys Biosciences
365 N 600 W
Logan, UT 84321
(435) 752-0531

Q-Plex™ Human Custom Singleplex

REF	Part Number: 107749GR-01	Software Product Code: HCUM250922 
LOT	Lot Number: HCUK250922 22700 001	
	Expiration Date: 2026-JUL-24	

Reconstitute Calibrator using 500 µL of Complete Sample Diluent

Component	Part Number	Lot Number	Recon. Volume	Exp. Date
Plate	101852HU	HCUP250925	N/A	2027-SEP-25
Sample Diluent	101122HU	HQ16S250825	N/A	2026-AUG-25
Competitor	101890GR	HCUC250912	N/A	2026-SEP-12
Calibrator	101855GR	HCUM250922	500 µL	2027-OCT-10
Detection	101856GR	HCUD250923	N/A	2026-SEP-23
Streptavidin HRP	101173GR	QNSN250124	N/A	2026-JUL-24
Substrate A	101159GR	QNSA250508	N/A	2026-DEC-09
Substrate B+	101120GR	QNSB+250508	N/A	2026-DEC-09
Wash Buffer	101158GR	QNSW250602	N/A	2027-JUN-02



Spot #	Analyte	Cal Conc	Units	LLD	LLOQ	ULOQ
1	PDG	2557	ng/mL	1.75	3.51	2557
2	REF	N/A	N/A	N/A	N/A	N/A

Q-Plex™ ARRAYS

QC Release:



Date:

13OCT25

QA Release:



Date:

13OCT25

for Research Use Only