

101896GR - 8.1

QUANSYS BIOSCIENCES
CERTIFICATE OF ANALYSIS


Quansys Biosciences
 365 N 600 W
 Logan, UT 84321
 (435) 752 - 0531

Q-Plex™ Human Cytokine Panel 1 (6-Plex)

REF	Part Number: 115233HU	Software Product Code: HA1M190516-P16P 
LOT	Lot Number: HCP1K190516 23058 001	
	Expiration Date: 2026-JUL-24	

Reconstitute Calibrator using 600 µL of Sample Diluent

Component	Part Number	Lot Number	Recon. Volume	Exp. Date
Plate	115252HU	HCP1P250508	N/A	2027-MAY-08
Sample Diluent	101104HU	HQ6S250205	N/A	2027-FEB-05
Calibrator	102655HU	HA1M190516	600 µL	2026-DEC-11
Detection	110956HU 01	HCSDQ6250818	N/A	2026-AUG-18
Streptavidin HRP	101173GR	QNSN250124	N/A	2026-JUL-24
Substrate A	101159GR	QNSA251013	N/A	2026-DEC-09
Substrate B+	101120GR	QNSB+251013	N/A	2026-DEC-09
Wash Buffer	101158GR	QNSW250128	N/A	2027-JAN-28



Spot #	Analyte	Cal Conc	Units	LLD	LLOQ	ULOQ
1	IL-1β	10759	pg/mL	7.38	14.76	10759
2	IL-6	2700	pg/mL	1.85	3.70	2700
3	TNFα	4912	pg/mL	3.37	6.74	4912
4	IFNγ	5700	pg/mL	3.91	7.82	5700
5	IL-1α	3400	pg/mL	2.33	4.66	3400
6	IL-10	6431	pg/mL	4.41	8.82	6431
7	REF	N/A	N/A	N/A	N/A	N/A

Q-Plex™ ARRAYS

QC Release:

Kaylee Kowalles

Date:

11 Dec 25

QA Release:

2 Cox

Date:

*12 DEC 25***for Research Use Only**