

101896GR - 8.1

QUANSYS
BIOSCIENCES
CERTIFICATE OF ANALYSIS


Quansys Biosciences
365 N 600 W
Logan, UT 84321
(435) 752-0531

Q-Plex™ Human Micronutrient v2 (7-Plex)

REF	Part Number: 355149HU	Software Product Code:
LOT	Lot Number: HMTK251113 23086 001	HMTM251113
	Expiration Date: 2026-SEP-11	

Reconstitute Calibrator using 270 µL of Complete Sample Diluent and 30 µL of Calibrator Additive

Component	Part Number	Lot Number	Recon. Volume	Exp. Date
Plate	355152HU	HMTM251203	N/A	2027-JUN-03
Sample Diluent	355157HU	HMTS251121	N/A	2027-MAY-21
Competitor	355190HU	HMTC251124	20 mL	2027-JUN-19
Calibrator	355155HU	HMTM251113	270 µL	2027-JUN-19
Calibrator Additive	355142HU	HMTA250311	N/A	2026-SEP-11
Detection	355156HU	HMTD251205	N/A	2027-JUN-05
Streptavidin HRP	101173GR	QNSN250610	N/A	2026-DEC-10
Substrate A	101159GR	QNSA251013	N/A	2026-DEC-09
Substrate B2	101128GR	QNSB2251013	N/A	2026-DEC-09
Wash Buffer	101158GR	QNSW251210	N/A	2027-DEC-10



Spot #	Analyte	Cal Conc	Units	LLD	LLOQ	ULOQ
1	Ferritin	22.60	µg/L	0.016	0.031	22.60
2	*AGP	0.36	g/L	0.00074	0.0015	0.36
3	Tg	14.34	µg/L	0.012	0.024	14.34
4	sTfR	3.51	mg/L	0.0024	0.0048	3.51
5	REF	N/A	N/A	N/A	N/A	N/A
6	CRP	0.96	mg/L	0.00066	0.0013	0.17
7	HRP2	1.44	µg/L	0.0051	0.0100	1.44
8	*RBP4	0.94	µmol/L	0.0019	0.0039	0.94

Q-Plex™ ARRAYS
*Competitive Assay

QC Release: Kayla Korallis Date: 22 Dec 25
 QA Release: [Signature] Date: 22 DEC 25
 for Research Use Only