



CERTIFICATE OF ANALYSIS


Q-Plex™ Human Cytokine HS (4-Plex)

REF	Part Number: 112533HU	Software Product Code: HCHM200803-4P 
LOT	Lot Number: HCH4K200803 23381 001	
	Expiration Date: 2026-JUL-24	

Reconstitute Calibrator using 750 µL of Sample Diluent

Component	Part Number	Lot Number	Recon. Volume	Exp. Date
Plate	112552HU	HCH4P260212	N/A	2028-FEB-12
Sample Diluent	101104HU	HQ6S250205	N/A	2027-FEB-05
Calibrator	112455HU	HCHM200803	750 µL	2028-FEB-25
Detection	112456HU 01	HCHDQ6260220	6 mL	2028-FEB-20
Streptavidin HRP	101173GR	QNSN250124	N/A	2026-JUL-24
Substrate A	101159GR	QNSA251218	N/A	2026-DEC-09
Substrate B+	101120GR	QNSB+251218	N/A	2026-DEC-09
Wash Buffer	101158GR	QNSW250425	N/A	2027-APR-25

Spot #	Analyte	Cal Conc	Units	LLD	LLOQ	ULOQ
1	IL-4	18.00	pg/mL	0.012	0.025	18.00
2	IL-6	240	pg/mL	0.16	0.33	240
3	IL-10	2320	pg/mL	1.59	3.18	2320
4	IFN γ	93.00	pg/mL	0.064	0.13	93.00
5	REF	N/A	N/A	N/A	N/A	N/A



QC Release: Kayla Kowallis Date: 25 FEB 26

QA Release: RL Cox Date: 25 FEB 26