



101896GR - 8.1

QUANSYS
BIOSCIENCES
CERTIFICATE OF ANALYSIS


Quansys Biosciences
365 N 600 W
Logan, UT 84321
(435) 752 - 0531

Q-Plex™ Human Custom 15-Plex

REF	Part Number: 107749GR-15	Software Product Code: HA1M251204-L15 DV 
LOT	Lot Number: HCUK251204 23238 001	
	Expiration Date: 2026-JUL-24	

Reconstitute Calibrator 1 using 650 µL of Sample Diluent
Reconstitute Calibrator 2 using 400 µL of Calibrator 1

Component	Part Number	Lot Number	Recon. Volume	Exp. Date
Plate	101852HU	HCUP260225	N/A	2028-FEB-25
Sample Diluent	101104HU	HQ6S250205	N/A	2027-FEB-05
Calibrator 1	102655HU	HA1M251204	650 µL	2028-MAR-09
Calibrator 2	102855HU	HA2M200303	400 µL	2028-MAR-09
Detection	101856HU	HCUD260302	N/A	2028-MAR-02
Streptavidin HRP	101173GR	QNSN250124	N/A	2026-JUL-24
Substrate A	101159GR	QNSA251218	N/A	2026-DEC-09
Substrate B+	101120GR	QNSB+251218	N/A	2026-DEC-09
Wash Buffer	101158GR	QNSW250425	N/A	2027-APR-25



Q-Plex™ ARRAYS

Spot #	Analyte	Cal Conc	Units	LLD	LLOQ	ULOQ
1	INFγ	8930	pg/mL	6.12	12.25	8930
2	IL-1α	4540	pg/mL	3.11	6.23	4540
3	IL-1β	13480	pg/mL	9.25	18.49	13480
4	IL-6	3210	pg/mL	2.20	4.40	3210
5	IL-8	2280	pg/mL	1.56	3.13	2280
6	IL-10	9050	pg/mL	6.21	12.41	9050
7	IL-12p70	3000	pg/mL	2.06	4.12	3000
8	IL-21	280	pg/mL	0.19	0.38	280
9	IP-10	2730	pg/mL	1.87	3.74	910
10	MCP-1	2180	pg/mL	4.49	8.97	727
11	MCP-2	590	pg/mL	3.64	7.28	590
12	MIP-1α	2400	pg/mL	1.65	3.29	800
13	MIP-1β	1140	pg/mL	2.35	4.69	1140
14	RANTES	6730	pg/mL	4.62	9.23	6730
15	TNFα	7100	pg/mL	4.87	9.74	7100
16	REF	N/A	N/A	N/A	N/A	N/A

QC Release: *Kaylee Kowalik* Date: 09MAR26
 QA Release: *[Signature]* Date: 09MAR26
for Research Use Only