



Quansys Biosciences, Inc. | 365 N 600 W Logan, UT 84321 | (435) 752-0531

CERTIFICATE OF ANALYSIS


Q-Plex™ Human Cytokine Release Syndrome (16-Plex)

REF	Part Number: 119733HU	Software Product Code: HCRM200421-CR2 
LOT	Lot Number: HCRK200421 23047 001	
EXP	Expiration Date: 2027-OCT-21	

Reconstitute Calibrator using 550µL of Sample Diluent

Component	Part Number	Lot Number	Recon. Volume	Exp. Date
Plate	119752HU	HCRP251217	N/A	2027-DEC-17
Sample Diluent	101104HU	HQ6S251021	N/A	2027-OCT-21
Calibrator	119755HU	HCRM200421	550µL	2028-APR-05
Detection	119756HU	HCRD260414	N/A	2028-APR-14
Streptavidin HRP	101173GR	QNSN260507	N/A	2027-NOV-07
Substrate A	101159GR	QNSA260427	N/A	2028-APR-02
Substrate B+	101120GR	QNSB+260427	N/A	2028-APR-02
Wash Buffer	101158GR	QNSW260227	N/A	2028-FEB-27

Spot #	Analyte	Cal Conc	Units	LLD	LLOQ	ULOQ
1	GMCSF	6900	pg/mL	4.73	9.47	6900
2	IFNα	390.00	pg/mL	0.27	0.53	390.00
3	IFNγ	95.00	pg/mL	0.070	0.13	95.00
4	IL-1β	3900	pg/mL	2.67	5.35	3900
5	IL-1Ra	7800	pg/mL	5.35	10.70	7800
6	IL-2	67.00	pg/mL	0.050	0.090	67.00
7	IL-2Ra	19500	pg/mL	13.37	26.75	19500
8	IL-6	360.00	pg/mL	0.25	0.49	360.00
9	IL-6R	1400	pg/mL	0.96	1.92	466.67
10	IL-8	2200	pg/mL	1.51	3.02	2200
11	IL-10	2500	pg/mL	1.71	3.43	2500
12	IL-12p70	3000	pg/mL	2.06	4.12	3000
13	IL-13	340.00	pg/mL	0.23	0.47	340.00
14	MCP-1	2600	pg/mL	5.35	10.70	866.67
15	MIP-1α	2200	pg/mL	1.51	3.02	733.33
16	TNFα	630.00	pg/mL	0.43	0.86	630.00
17	REF	N/A	N/A	N/A	N/A	N/A



QC Release: Kymle Kowallis Date: 08 JUN 26

QA Release: [Signature] Date: 08 JUN 26

101896GR - 9.0

For Research Use Only