

CERTIFICATE OF ANALYSIS


Q-Plex™ Human Micronutrient v2 (7-Plex)

REF	Part Number: 355149HU	Software Product Code: HMTMA260430 23742 001 
LOT	Lot Number: HMTK251113 23742 001	
	Expiration Date: 2027-OCT-30	

Reconstitute Calibrator using 270 µL of Complete Sample Diluent and 30 µL of Calibrator Additive

Component	Part Number	Lot Number	Recon. Volume	Exp. Date
Plate	355152HU	HMTTP260505	N/A	2027-NOV-05
Sample Diluent	355157HU	HMTS260505	N/A	2027-NOV-05
Competitor	355190HU	HMTC251124	20 mL	2027-DEC-03
Calibrator	355155HU	HMTM251113	270 µL	2027-DEC-03
Calibrator Additive	355142HU	HMTA260430	N/A	2027-OCT-30
Detection	355156HU	HMTD260514	N/A	2027-NOV-14
Streptavidin HRP	101173GR	QNSN260507	N/A	2027-NOV-07
Substrate A	101159GR	QNSA260427	N/A	2028-APR-02
Substrate B2	101128GR	QNSB2260427	N/A	2028-APR-02
Wash Buffer	101158GR	QNSW260528	N/A	2028-MAY-28

Spot #	Analyte	Cal Conc	Units	LLD	LLOQ	ULOQ
1	Ferritin	26.66	µg/L	0.020	0.039	16.52
2	*AGP	0.36	g/L	0.00031	0.00063	0.32
3	Tg	14.34	µg/L	0.012	0.025	11.46
4	sTfR	3.51	mg/L	0.0035	0.0070	2.84
5	REF	N/A	N/A	N/A	N/A	N/A
6	CRP	0.96	mg/L	0.00064	0.0013	0.17
7	HRP2	1.44	µg/L	0.0051	0.0100	1.44
8	*RBP4	0.94	µmol/L	0.00079	0.0016	0.72



*Competitive Assay

QC Release:

Kaylee Howells

Date:

04 JUN 26

QA Release:

[Signature]

Date:

04 JUN 26

101896GR - 9.0

for Research Use Only