



CERTIFICATE OF ANALYSIS


Q-Plex™ Human Pneumococcal IgG (9-Plex) 30 mL Diluent

REF	Part Number: 482149HU-30	Software Product Code: HPNM260514 23763 001 
LOT	Lot Number: HPNK260514 23763 001	
	Expiration Date: 2027-MAY-08	

Reconstitute Calibrator using 500µL of Sample Diluent

Component	Part Number	Lot Number	Recon. Volume	Exp. Date
Plate	482152HU	HPNP260527	N/A	2027-MAY-27
Sample Diluent	481157HU-30	HPNS260529	N/A	2027-MAY-29
Calibrator	481155HU	HPNM260514	500µL	2027-JUN-02
Detection	482156HU	HPND260508	N/A	2027-MAY-08
Substrate A	101159GR	QNSA260427	N/A	2028-APR-02
Substrate B+	101120GR	QNSB+260427	N/A	2028-APR-02
Wash Buffer	101158GR	QNSW260518	N/A	2028-MAY-18

Spot #	Analyte	Cal Conc	Units	LLD	LLOQ	ULOQ
1	CWPS	N/A	N/A	N/A	N/A	N/A
2	2	2.45	µg/mL	0.00063	0.0013	1.23
3	10A	1.28	µg/mL	0.00044	0.00088	0.45
4	11A	0.46	µg/mL	0.00014	0.00029	0.14
5	15B	1.63	µg/mL	0.00058	0.0012	0.60
6	17F	0.85	µg/mL	0.00044	0.00089	0.52
7	19A	1.42	µg/mL	0.00065	0.0013	1.31
8	20	0.99	µg/mL	0.00034	0.00068	0.59
9	22F	1.06	µg/mL	0.00042	0.00083	0.64
10	33F	1.08	µg/mL	0.00052	0.0010	0.35
11	Control Spot	N/A	N/A	N/A	N/A	N/A
12	REF	N/A	N/A	N/A	N/A	N/A



QC Release: Kayla Kowallis Date: 02 JUN 26
 QA Release: [Signature] Date: 02 JUN 26