

CERTIFICATE OF ANALYSIS


Q-Plex™ Mouse Cytokine Panel 3 HS (4-Plex)

REF	Part Number: 111249MS	Software Product Code: MHSM230927-P34P 
LOT	Lot Number: MCP3K230927 23880 001	
	Expiration Date: 2026-NOV-07	

Reconstitute Calibrator using 550µL of Sample Diluent

Component	Part Number	Lot Number	Recon. Volume	Exp. Date
Plate	111252MS	MCP3P260529	N/A	2027-MAY-29
Sample Diluent	112457MS	MHSS250917	N/A	2027-MAY-17
Calibrator	112455MS	MHSM230927	550µL	2027-JUN-10
Detection	112456MS	MHSD251107	6 mL	2026-NOV-07
Streptavidin HRP	101173GR	QNSN251216	N/A	2027-JUN-16
Substrate A	101159GR	QNSA260521	N/A	2026-DEC-09
Substrate B+	101120GR	QNSB+260521	N/A	2026-DEC-09
Wash Buffer	101158GR	QNSW260528	N/A	2028-MAY-28

Spot #	Analyte	Cal Conc	Units	LLD	LLOQ	ULOQ
1	IL-2	180	pg/mL	0.12	0.25	180
2	IL-4	160	pg/mL	0.11	0.22	160
3	IL-6	1110	pg/mL	0.76	1.52	1110
4	TNFα	390	pg/mL	0.27	0.53	390
5	REF	N/A	N/A	N/A	N/A	N/A



QC Release:

Kayla Kowallis

Date: 10JUN26

QA Release:

Myra Gal

Date: 10JUN26