



Quansys Biosciences, Inc. | 365 N 600 W Logan, UT 84321 | (435) 752-0531

CERTIFICATE OF ANALYSIS


Q-Plex™ Human MMP (6-Plex)

REF	Part Number: 340949HU	Software Product Code: MMPM220518
LOT	Lot Number: MMPK220518 23986 001	
EXP	Expiration Date: 2026-DEC-09	

Reconstitute Calibrator using 300 µL of Complete Sample Diluent

Component	Part Number	Lot Number	Recon. Volume	Exp. Date
Plate	340952HU	MMPP250721	N/A	2027-JUL-21
Sample Diluent	340957HU	MMPS250721	N/A	2027-JUL-21
Assay Diluent	340983HU	MMPV250722	N/A	2027-JUL-22
Calibrator	340955HU	MMPM220518	300 µL	2028-JUN-10
Detection	340956HU	MMPD250729	N/A	2027-JUL-29
Streptavidin HRP	101173GR	QNSN251216	N/A	2027-JUN-16
Substrate A	101159GR	QNSA260226	N/A	2026-DEC-09
Substrate B+	101120GR	QNSB+260226	N/A	2026-DEC-09
Wash Buffer	101158GR	QNSW260411	N/A	2028-APR-11

Spot #	Analyte	Cal Conc	Units	LLD	LLOQ	ULOQ
1	MMP 1	12563	pg/mL	25.85	51.70	12563
2	MMP 2	85197	pg/mL	175.30	350.60	85197
3	MMP 3	4579	pg/mL	9.42	18.84	4579
4	MMP 7	1745	pg/mL	3.59	7.18	1745
5	MMP 9	44914	pg/mL	92.42	184.83	44914
6	MMP 13	32943	pg/mL	67.78	135.57	32943
7	REF	N/A	N/A	N/A	N/A	N/A



QC Release:

K. Cox

Date:

10 JUN 26

QA Release:

Kaylee Howells

Date:

10 JUN 26

101896GR - 9.0

for Research Use Only