



Quansys Biosciences, Inc. | 365 N 600 W Logan, UT 84321 | (435) 752-0531

CERTIFICATE OF ANALYSIS


Q-Plex™ Human Pneumococcal IgG Custom 13-Plex

REF	Part Number: 107749GR-13	Software Product Code: HPNM240514-O13 
LOT	Lot Number: HCUK240514 24042 001	
EXP	Expiration Date: 2026-DEC-09	

Reconstitute Calibrator using 500µL of Sample Diluent

Component	Part Number	Lot Number	Recon. Volume	Exp. Date
Plate	101852GR	HCUP260629	N/A	2027-JUN-29
Sample Diluent	481157HU	HPNS260624	N/A	2027-JUN-24
Calibrator	481155HU	HPNM240514	500µL	2027-JUL-09
Detection	482156HU	HPND260119	N/A	2027-JAN-19
Substrate A	101159GR	QNSA260521	N/A	2026-DEC-09
Substrate B+	101120GR	QNSB+260521	N/A	2026-DEC-09
Wash Buffer	101158GR	QNSW260707	N/A	2028-JUL-07

Spot #	Analyte	Cal Conc	Units	LLD	LLOQ	ULOQ
1	1	0.76	µg/mL	0.00037	0.00074	0.19
2	3	0.12	µg/mL	0.00023	0.00047	0.12
3	4	0.27	µg/mL	0.00013	0.00026	0.068
4	5	0.65	µg/mL	0.00032	0.00063	0.16
5	6B	0.73	µg/mL	0.00036	0.00071	0.18
6	7F	0.73	µg/mL	0.0014	0.0029	0.18
7	8	1.27	µg/mL	0.0025	0.0050	0.32
8	9N	0.59	µg/mL	0.0012	0.0023	0.15
9	9V	0.58	µg/mL	0.0011	0.0023	0.15
10	14	3.49	µg/mL	0.00039	0.00078	0.050
11	18C	0.65	µg/mL	0.0017	0.0034	0.22
12	19F	1.27	µg/mL	0.0013	0.0025	0.16
13	23F	0.57	µg/mL	0.0025	0.0050	0.32
14	CWPS	N/A	N/A	N/A	N/A	N/A
15	Control Spot	N/A	N/A	N/A	N/A	N/A
16	REF	N/A	N/A	N/A	N/A	N/A



Q-Plex™ ARRAYS

QC Release:

Kayla Kowallis

Date: 10 JUL 26

QA Release:

Na Ja

Date: 10 JUL 26

101896GR - 9.0

For Research Use Only