

CERTIFICATE OF ANALYSIS

Q-Plex™ Human Micronutrient v2 (7-Plex)


REF	Part Number: 355149HU	Software Product Code:
LOT	Lot Number: HMTK251113 24078 001	HMTMA260430 24078 001
	Expiration Date: 2027-JUL-09	

Reconstitute Calibrator using 270 µL of Complete Sample Diluent and 30 µL of Calibrator Additive

Component	Part Number	Lot Number	Recon. Volume	Exp. Date
Plate	355152HU	HMTTP260630	N/A	2027-DEC-30
Sample Diluent	355157HU	HMTS260109	N/A	2027-JUL-09
Competitor	355190HU	HMTC251124	20 mL	2028-JAN-13
Calibrator	355155HU	HMTM251113	270 µL	2028-JAN-13
Calibrator Additive	355142HU	HMTA260430	N/A	2027-OCT-30
Detection	355156HU	HMTD260514	N/A	2027-NOV-14
Streptavidin HRP	101173GR	QNSN260507	N/A	2027-NOV-07
Substrate A	101159GR	QNSA260612	N/A	2028-APR-02
Substrate B2	101128GR	QNSB2260612	N/A	2028-APR-02
Wash Buffer	101158GR	QNSW260604	N/A	2028-JUN-04

Spot #	Analyte	Cal Conc	Units	LLD	LLOQ	ULOQ
1	Ferritin	26.66	µg/L	0.020	0.039	16.52
2	*AGP	0.36	g/L	0.00031	0.00063	0.32
3	Tg	14.34	µg/L	0.012	0.025	11.46
4	sTfR	3.51	mg/L	0.0035	0.0070	2.84
5	REF	N/A	N/A	N/A	N/A	N/A
6	CRP	0.96	mg/L	0.00064	0.0013	0.17
7	HRP2	1.44	µg/L	0.0051	0.010	1.44
8	*RBP4	0.94	µmol/L	0.00079	0.0016	0.72

*Competitive Assay



Q-Plex™ ARRAYS

QC Release:

Kayla Konvalis

Date: 13 JUL 26

QA Release:

[Signature]

Date: 13 JUL 26

101896GR - 9.0

For Research Use Only