

CERTIFICATE OF ANALYSIS


Q-Plex™ Human Malaria (5-Plex)

REF	Part Number: 565949HU	Software Product Code: HMRM191125-ML2 
LOT	Lot Number: HMRK191125 24251 001	
	Expiration Date: 2027-JUL-01	

Reconstitute Calibrator using 350µL of Complete Sample Diluent

Component	Part Number	Lot Number	Recon. Volume	Exp. Date
Plate	565952HU	HMRP260601	N/A	2027-DEC-01
Calibrator Diluent	565928HU	HMRU260714	6mL	2028-JAN-14
Assay Diluent	565983HU	HMRV260618	N/A	2027-DEC-18
Competitor	565990HU	HMRC191014	6mL	2028-FEB-14
Calibrator	565955HU	HMRM191125	350µL	2028-FEB-14
Detection	565956HU	HMRD260715	N/A	2028-JAN-15
Streptavidin HRP	101173GR	QNSN260507	N/A	2027-NOV-07
Substrate A+	105959GR	MOSA+260701	N/A	2027-JUL-01
Substrate B	105920GR	MOSB260701	N/A	2027-JUL-01
Wash Buffer	101158GR	QNSW260518	N/A	2028-MAY-18

Spot #	Analyte	Cal Conc	Units	LLD	LLOQ	ULOQ
1	HRP2	2800	pg/mL	0.53	0.68	2800
2	LDH-Pan	67000	pg/mL	7.28	16.36	67000
3	LDH-Pv	19200	pg/mL	4.23	4.69	19200
4	LDH-Pf	20800	pg/mL	2.33	5.08	20800
5	*CRP	38000	ng/mL	7.65	9.28	38000
6	REF	N/A	N/A	N/A	N/A	N/A



*Competitive Assay

QC Release:

Country Namb

Date:

14Aug26

QA Release:

H. Lax

Date:

14Aug26

101896GR - 9.0

for Research Use Only